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Material Safety Data Sheet (MSDS) for passengerships and airlines.

Shenzhen SEM.Test Technology Co., Ltd.  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

Material Safety Data Sheet

Material Safety Data Sheet (MSDS)

For  
Changzhou Hankun New Energy Technology Co., Ltd.  
No. 2188, Longcheng Avenue, Zhonglou District, Changzhou city, Jiangsu province, China  
and for their product

Polymer lithium-ion battery pack

Model/type reference ..... HKE12-10  
Trademark ..... HanKun  
Nominal Voltage..... 12.0V  
Typical Capacity..... 10.0Ah, 120.0Wh  
Weight..... 1294.0g  
Shape and Physical Dimension L: 505.0mm  
(mm)..... W: 68.6mm  
T: 31.6mm  
Version number..... V3.0  
Preparation Date..... Feb. 05 , 2018  
Revision date..... N/A.

Laboratory ..... Shenzhen SEM.Test Technology Co., Ltd.  
Address ..... 1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

Compiled by (name+ signature) .. Ivan Ni

Ivan Ni

Approved by (name+ signature) .. Davis Yang

Davis yang

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

### Product Identifier

Product name: Polymer lithium-ion battery pack

Model: HKE12-10

### Other means of identification

Synonyms:none

### Recommended use of the chemical and restrictions on use

Recommended Use:Used in portable electronic equipments;

Uses advised against:

- a) Do not dismantle, open or shred secondary cells or batteries.
- b) Do not expose cells or batteries to heat or fire. Avoid storage in direct sunlight.
- c) Do not short-circuit a cell or a battery. Do not store cells or batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- d) Do not remove a cell or battery from its original packaging until required for use.
- e) Do not subject cells or batteries to mechanical shock.
- f) In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- g) Do not use any charger other than that specifically provided for use with the equipment.
- h) Observe the plus (+) and minus (-) marks on the cell, battery and equipment and ensure correct use.
- i) Do not use any cell or battery which is not designed for use with the equipment.
- j) Do not mix cells of different manufacturer, capacity, size or type within a device.
- k) Battery usage by children should be supervised.
- l) Seek medical advice immediately if a cell or a battery has been swallowed.
- m) Always purchase the battery recommended by the device manufacturer for the equipment.
- n) Keep cells and batteries clean and dry.
- o) Wipe the cell or battery terminals with a clean dry cloth if they become dirty.
- p) Secondary cells and batteries need to be charged before use. Always use the correct charger and refer to the manufacturer's instructions or equipment manual for proper charging instructions.
- q) Do not leave a battery on prolonged charge when not in use.
- r) After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
- s) Retain the original product literature for future reference.
- t) Use only the cell or battery in the application for which it was intended.
- u) When possible, remove the battery from the equipment when not in use.
- v) Dispose of properly.

### Details of the supplier of the safety data sheet:

Supplier Name: Changzhou Hankun New Energy Technology Co., Ltd.

Address: No. 2188, Longcheng Avenue, Zhonglou District, Changzhou city, Jiangsu province, China

Telephone number of the supplier: 0086-0519-83383316

Emergency Telephone No.(24h): 0086-0519-83383316

Fax: 0086-0519-83383316

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

Postcode:213012

E-mail address: shane.zhao@163.com

**Emergency telephone number**

**Company Emergency Phone Number:** 0086-0519-83383316

## 2. HAZARDS IDENTIFICATION

### Classification

No harm at the normal use. If contact the Electrolyte liquid in the Lithium ion battery, reference as follows:

### Classification of the substance or mixture

Classification according to GHS

Acute Toxicity, Oral(Hazard category 4)

Acute Toxicity, Dermal(Hazard category 3)

Skin, irritate(Cagegory 1B)

Eye Irritate (Hazard category 1)

### GHS Label elements, including precautionary statements:

GHS02

GHS05

GHS06

**Signal word:** Warning

**Hazard statement(s):**

H242:Heating may cause a fire;

H311: Toxic in contact with skin;

H314:Causes severe skin burns and eye damage;

H302:Harmful if swallowed;

### precautionary statements:

**Prevention:**

P264 Wash thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

**Response:**

P312:Call a Poison center or doctor/physician if you feel unwell.

P302+P350-IF ON SKIN: Gently wash with plenty of soap and water

P301+P330+P331-IF SWALLOWED: rise mouth. Do NOT induce vomiting

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**Storage:**

None

**Disposal**

P501: Dispose of contents/container in accordance with local/national regulations

**Hazards not otherwise classified (HNOC)**

Not Applicable

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

**Other information**

No information available.

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

**Chemical characterization:** Mixtures

**Description:**

Product: Consisting of the following components.

Common Chemical Name	Concentration (%)	CAS Number	EC No.
Aluminum Foil (Al)	5	7429-90-5	231-072-3
Poly Vinylidene Fluoride PVDF(-[CH <sub>2</sub> -CF <sub>2</sub> ]-n)	6	24937-79-9	----
Graphite (C)	8	7782-42-5	231-955-3
Copper Foil (Cu)	7-10	7440-50-8	231-159-6
dimethyl carbonate (DMC)	2-5	616-38-6	----
Nickel	2-5	7440-02-0	----
Carbon(C)	12-15	7400-44-0	----
Ethylene carbonate	5	96-49-1	202-510-0
LiFePO <sub>4</sub>	35-40	N/A	----
Li	2-3	7439-93-2	----

Note: CAS number is Chemical Abstract Service Registry Number.

N/A=Not apply.

### 4. FIRST-AID MEASURES

**First aid measures**

**Eye Contact** Rinse thoroughly with plenty of water, also under the eyelids. If symptoms persist, call a physician.

**Skin Contact** Remove contaminated clothing and shoes. Wash skin with soap and water. In the case of skin irritation or allergic reactions see a physician.

**Inhalation** Move to fresh air. If symptoms persist, call a physician.

**Ingestion** Do NOT induce vomiting. Drink plenty of water. If symptoms persist, call a physician.

Most important symptoms and effects, both acute and delayed

Swallowing Do not induce vomiting. Get medical attention.

**Most Important Symptoms/Effects** No information available.

**Indication of any immediate medical attention and special treatment needed**

Notes to Physician Treat symptomatically

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

*Material Safety Data Sheet*

## 5. FIRE-FIGHTING MEASURES

### Suitable Extinguishing Media

CO<sub>2</sub>, dry chemical powder, water spray.

Unsuitable Extinguishing Media: No information available.

### Specific Hazards Arising from the Chemical

Formation of toxic gases is possible during heating or in case of fire.

In case of fire, the following can be released:

Carbon monoxide(CO)

Carbon dioxide

Other irritating and toxic gases.

### Hazardous Combustion Products

Carbon oxides.

Explosion Data

Sensitivity to Mechanical Impact No

Sensitivity to Static Discharge No

### Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear. For example: Wear self-contained respiratory protective device. Wear suitable protective clothing and eye/face protection.

### Special hazards arising from the substance or mixture:

Battery may burst and release hazardous decomposition products when exposed to a fire situation.

Lithium ion batteries contain flammable electrolyte that may vent, ignite and produce sparks when subjected to high temperature(>150°C). When damaged or abused(e.g. mechanical damage or electrical overcharging); may burn rapidly with flare-burning effect; may ignite other batteries in clothes proximity.

## 6. ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

Personal Precautions Avoid contact with eyes.

Refer to section 8 for personal protective equipment. Ensure adequate ventilation. Remove all sources of ignition.

Evacuate personnel to safe areas.

### Environmental precautions

Environmental Precautions Refer to protective measures listed in Sections 7 and 8.

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

### Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Cleaning up Use personal protective equipment. Dam up. Cover liquid spill with sand, earth or other Non combustible absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly.

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

## 7. HANDLING AND STORAGE

### Precautions for safe handling

Handling Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes and clothing. Wear personal protective equipment.  
Wash thoroughly after handling. Use this material with adequate ventilation.  
The product is not explosive.

### Conditions for safe storage, including any incompatibilities

If the Polymer lithium-ion battery pack is subject to storage for such a long term as more than 3 months, it is recommended to recharge the Polymer lithium-ion battery pack periodically.  
3 months: -10°C~+40°C, 45 to 85%RH

And recommended at 0°C~+35°C for long period storage.

The capacity recovery rate in the delivery state (50% capacity of fully charged) after storage is assumed to be 80% or more.

The voltage for a long time storage shall be 11.0V~12.8V range.

Do not storage Polymer lithium-ion battery pack haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.

Keep out of reach of children.

Do not expose Polymer lithium-ion battery pack to heat or fire. Avoid storage in direct sunlight.

Do not store together with oxidizing and acidic materials.

Keep ignition sources away- Do not smoke.

Store in cool, dry and well-ventilated place.

Incompatible Products None known.

## 8. EXPOSURE CONTROLS AND PERSONAL PROTECTION

### Control parameters

Ingredients with limit values that require monitoring at the workplace:

12190-79-3 Lithium Cobalt Oxide

TLV (USA)	0.02mg/m <sup>3</sup>
MAK (Germany)	0.1mg/m <sup>3</sup>

Other Exposure Guidelines Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962(11th Cir., 1992).

### Appropriate engineering controls

Engineering Measures Showers

Eyewash stations

Ventilation systems

Use adequate general or local exhaust ventilation to keep airborne concentrations below the permissible exposure limits. Ensure adequate ventilation.

### Individual protection measures, such as personal protective equipment

Eye/Face Protection:

Shenzhen SEM.Test Technology Co., Ltd.

1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

Material Safety Data Sheet

**Tightly sealed goggles**

**Body protection:**

Protective work clothing.

**Skin protection:**

**Protective gloves**

**Material of gloves:**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

**Penetration time of glove material:**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

**Respiratory Protection** No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Physical State	Form: Prismatic
	Color: Black
	Odour: Odourless
	Odor Threshold: No information available
Change in condition:	
pH, with indication of the concentration	Not determined.
Melting point/freezing point	Not determined.
Initial boiling point and Boiling range:	Not determined.
Flash Point	Not determined.
Evaporation rate	Not determined.
Flammability (solid, gas)	Not determined.
Upper/lower flammability or explosive limits	Not determined.
Vapor Pressure:	Not determined.
Vapor Density:	Not determined.
relative density:	Not determined.

Shenzhen SEM.Test Technology Co., Ltd.

1/F, Building A, Hongwei Industrial Park, Luxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

Material Safety Data Sheet

Solubility in Water:	Not determined.
Solubility in other solvents	Not determined.
n-octanol/water partition coefficient	Not determined.
Auto-ignition temperature	Product is not self-igniting.
Decomposition temperature	Not determined.
Odour threshold	Not determined.
Evaporation rate	Not determined.
Viscosity	Not determined.
Other Information	No further relevant information available.

## 10. STABILITY AND REACTIVITY

**Reactivity:** Stable under recommended storage and handling conditions (see section 7, Handling and storage).

**Chemical stability:** Stable under normal conditions of use, storage and transport.

**Thermal decomposition/conditions to be avoided:** No decomposition if used according to specifications.

**Possibility of Hazardous Reactions:** None under normal processing.

**Hazardous Polymerization:** Hazardous polymerization does not occur.

**Conditions to avoid:** Strong heating, fire, Incompatible materials.

**Incompatible materials:** Strong oxidizing agents, Strong acids, Base metals.

**Hazardous Decomposition Products:** Carbon oxides, Other irritating and toxic gases.

## 11. TOXICOLOGICAL INFORMATION

**Acute toxicity:** No data available.

LD/LC50 values relevant for classification:

Not available.

**Skin corrosion/irritation:** No irritant effect.

**Serious eye damage/irritation:** Cause serious eye irritation.

**Respiratory or skin sensitization:** No sensitizing effects known.

**Specific target organ system toxicity:** No information available.

**CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction):** No information available.

## 12. Ecological Information

**Toxicity:**

Acquatic toxicity:

**Shenzhen SEM.Test Technology Co., Ltd.**  
1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

No further relevant information available.

**Persistence and degradability:** No further relevant information available.

**Bioaccumulative potential:** No further relevant information available.

**Mobility in soil:** No further relevant information available.

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

**Other adverse effects:** No information available.

### 13. DISPOSAL CONSIDERATIONS

#### **Waste treatment methods**

Recommendation: Must not be disposed together with household garbage.

Do not allow product to reach sewage system

#### **Uncleaned packaging:**

Recommendation: Disposal must be made according to official regulations.

### 14. TRANSPORT INFORMATION

The Polymer lithium-ion battery pack had been tested according to the requirements of the UN manual of tests and Criteria, Part III, subsection 38.3;

The Polymer lithium-ion battery pack with a Watt-hour rating exceeding 100Wh or the cell with a Watt-hour rating in exceeding of 20Wh, The lithium ion batteries according to Section IA of PACKING INSTRUCTION 965, or Section I of PACKING INSTRUCTION 966~967 of the 2018 Dangerous Goods regulations 59th Edition may be transported.

The packaging shall be adequate to avoid mechanical damage during transport, handling and stacking. The materials and pack design shall be chosen so as to prevent the development of unintentional electrical conduction, corrosion of the terminals and ingress of moisture.

According to the Packing Instruction of IATA DGR 59th Edition for transportation.

Meets requirements of International Maritime Dangerous Goods(IMDG)-2016 Special Provision 188 to be transported as dangerous goods;

Meets the requirements of 49CFR173.185 to be transported as dangerous goods for road, rail, air, and vessel.

Meets the requirements of TDG special provision 34 to be transported as dangerous goods.

The package must be handled with care and that a flammability hazard exists if the package is damaged;

UN number of lithium battery: UN3480 or UN3481;

UN Proper shipping name/Description (technical name): Lithium ion batteries or Lithium ion batteries contained in equipment or Lithium ion batteries packed with equipment;

UN Classification (Transport hazard class): Non dangerous;

### 15. REGULATORY INFORMATION

#### **Safety, health and environmental regulations/legislation specific for the substance or mixture**

Report No.:SRT18029012S

Page 9 of 10

**Shenzhen SEM.Test Technology Co., Ltd.**  
 1/F, Building A, Hongwei Industrial Park, Liuxian 2nd  
 Road, Bao'an District, Shenzhen, P.R.C. (518101)

**Material Safety Data Sheet**

EU Regulation:

**Authorisations:** No information available.

**Restrictions on use:** No information available.

**Regulatory information**

CAS No.	EU (EINECS)	US (TSCA)	Japan (ENCS)	Canada (DSL/ NDSL)	Australia (AICS)	Korea (ECL)	China (IECSC)
7429-90-5	Listed	Not listed	Not listed	NDSL	Not listed	Not listed	Not listed
24937-79-9	Listed	Listed	Listed	DSL	Listed	Listed	Listed
7782-42-5	Not listed	Listed	Not listed	DSL	Listed	Listed	Listed
7440-50-8	Not listed	Listed	Not listed	DSL	Listed	Listed	Listed
616-38-6	Listed	Not listed	Not listed	NDSL	Not listed	Not listed	Not listed
7440-02-0	Listed	Listed	Listed	DSL	Listed	Listed	Listed
7400-44-0	Not listed	Listed	Not listed	DSL	Listed	Listed	Listed
96-49-1	Not listed	Listed	Not listed	DSL	Listed	Listed	Listed
7439-93-2	Listed	Not listed	Not listed	NDSL	Not listed	Not listed	Not listed

**Chemical safety assessment** A Chemical Safety Assessment has not been carried out.

## 16. OTHER INFORMATION

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

**Relevant phrases:**

R20/22: Harmful by inhalation and if swallowed.

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Name of Sample: Polymer Li-ion Battery HKE12-10 12V 10Ah 120Wh

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Shanghai Research Institute of Chemical Industry Testing Co., Ltd.

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Shanghai Research Institute of Chemical Industry  
Testing Co., Ltd. Test Report

NO. 1120010082

1/11

样品名称 Name of Sample	中文 Chinese	聚合物锂离子电池组 HKE12-10 12V 10Ah 120Wh	
	英文 English	Polymer Li-ion Battery HKE12-10 12V 10Ah 120Wh	
样品编号 Sample No.	1120010082		
委托单位 Consignor	常州瀚琨光伏新能源科技有限公司 Changzhou Hankun New Energy Technology Co., Ltd.		
生产单位 Manufacturer	常州瀚琨光伏新能源科技有限公司 Changzhou Hankun New Energy Technology Co., Ltd.		
检测方法 Test method	联合国《关于危险货物运输的建议书 试验和标准手册》 ST/SG/AC.10/11/Rev.6 Amend.1 38.3 UNITED NATIONS "Recommendations on the TRANSPORT OF DANGEROUS GOODS" Manual of Tests and Criteria ST/SG/AC.10/11/Rev.6 Amend.1 Section 38.3		
判定标准 Criterion	联合国《关于危险货物运输的建议书 试验和标准手册》 ST/SG/AC.10/11/Rev.6 Amend.1 38.3 UNITED NATIONS "Recommendations on the TRANSPORT OF DANGEROUS GOODS" Manual of Tests and Criteria ST/SG/AC.10/11/Rev.6 Amend.1 Section 38.3		
样品外观 Appearance	黑色 塑胶外壳 Black Plastics cement shell		
样品接受日期 Accepted Date	2020-01-06	检测起送日期 Test Date	2020-02-12 ~ 2020-02-28
	检测项目 Test Items		
检测项目 Test Items	高度模拟;热测试;振动;冲击;外短路;挤压;过充电;强制放电 Altitude simulation, Thermal test, Vibration, Shock, External short circuit, Crush, Overcharge, Forced discharge		
检测结论 Conclusion	经检测, 该样品符合联合国《关于危险货物运输的建议书 试验和标准手册》 ST/SG/AC.10/11/Rev.6 Amend.1 38.3标准要求。 The sample has passed the test items of UNITED NATIONS "Recommendations on the TRANSPORT OF DANGEROUS GOODS" Manual of Tests and Criteria ST/SG/AC.10/11/Rev.6 Amend.1 38.3		
备注 Comment	可充电锂电池组 Rechargeable Lithium Battery 生效日期(Date): 2020-03-02		
委托单位地址 Consignor Address	/		邮政编码 Post Code
			213002

批准  
Approver:  
职务  
Title:

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副总工程师(Vice chief engineer)

审核  
Checker:

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编制  
Compiler:

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# MSDS

## MATERIAL SAFETY DATA SHEET

**Prepared For** : Changzhou FengYu New Energy Technology Co., Ltd.  
No.35, Ieshan Road, Xinbei District, Changzhou City,  
Jiangsu Province, China

**Prepared By** : Shenzhen LCS Compliance Testing Laboratory Ltd.  
Room 101, 201, Building A and Room 301, Building C,  
Juji Industrial Park, Yabianxueziwei, Shajing Street,  
Bao'an District, Shenzhen, Guangdong, China

**Issue Date** : 2022.01.01

**Report Number** : KA2112161149AD100

**Written by:** Taylor Zhang

**Approved by:** \_\_\_\_\_

## Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

\* The MSDS is prepared based on the information provided by client. The contents and formats of this MSDS are revised as per client's request.

## Section 1-Chemical Product and Company Identification

<b>Product Name</b>	Polymer Lithium-Ion Battery
<b>Model</b>	FY24-8
<b>Trade Mark</b>	N/A
<b>Ratings</b>	24V, 8Ah,192Wh
<b>Weight</b>	1431.1g
<b>Manufacturer</b>	Changzhou FengYu New Energy Technology Co., Ltd.
<b>Manufacturer address</b>	No.35, Ieshan Road, Xinbei District, Changzhou City, Jiangsu Province, China
<b>Emergency Telephone</b>	+86-519-83859198
<b>Fax</b>	+86-519-83859198

## Section 2- Composition Information

Chemical Composition	CAS No.	Weight (%)	Trade Secret
Lithium cobaltate	12190-79-3	15 - 40	*
Graphite	7782-42-5	10 - 30	*
Phosphate(1-), hexafluoro-, lithium	21324-40-3	10 - 30	*
Copper	7440-50-8	7-13	*
Aluminium	7429-90-5	5-10	*
Nickel	7440-02-0	1-5	*

\*\*\* The exact percentage (concentration) of composition has been withheld as a trade secret.

## Section 3- Hazards Identification

<b>Emergency overview:</b>	N/A
<b>Classification according to GHS</b>	Not a dangerous substance according to GHS
<b>Label elements:</b>	
Hazard pictogram(s)	Not Applicable
Signal word	Not Applicable
Hazard statement(s)	Not Applicable

## Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

<b>Precautionary statement(s):</b>	
Prevention	Not Applicable
Response	Not Applicable
Disposal	Not Applicable
<b>Environmental hazards:</b>	No relevant information
<b>Important symptoms:</b>	See section 11 for more information
<b>Section 4- First Aid Measures</b>	
<b>Eye contact</b>	Flush eyes with plenty of water for least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical aid.
<b>Skin contact</b>	Remove contaminated clothes and rinse skin with plenty of water or shower for 15 minutes. Get medical aid.
<b>Inhalation</b>	Remove from exposure and move to fresh air immediately. Use oxygen if available.
<b>Ingestion</b>	Give at least 2 glasses of milk or water. Induce vomiting unless patient is unconscious. Call a physician.
<b>Section 5- Fire Fighting Measures</b>	
<b>Flash Point</b>	N/A
<b>Auto-Ignition Temperature</b>	N/A
<b>Extinguishing Media</b>	H <sub>2</sub> O, CO <sub>2</sub>
<b>Special Fire-Fighting Procedures</b>	Self-contained breathing apparatus
<b>Unusual Fire and Explosion Hazards</b>	Cell may vent when subjected to excessive heat-exposing battery contents
<b>Hazardous Combustion Products</b>	Carbon monoxide, carbon dioxide, lithium oxide fumes.
<b>Section 6- Accidental Release Measures</b>	
<b>Personal precautions, protective equipment and emergency procedures:</b> If the battery is released, remove personnel from area until fumes dissipate. Provide maximum ventilation to clear out hazardous gases. The preferred response is to leave the area and allow the vapors to dissipate. Avoid skin and eyes contact or inhalation of vapors. Remove spilled liquid with absorbent and incinerated. If leakage of the battery happens, liquid could be absorbed with sand, earth or other inert substance and contaminated area should be ventilated meantime.	
<b>Environment precautions:</b> Do not allow product to reach sewage system or any water source. Inform respective authorities in case of seepage into water course or sewage system. Do not allow to enter sewers/ surface or ground water.	

## Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

<b>Methods and material for containment and cleaning up:</b> If battery casing is dismantled, small amounts of electrolyte may leak. Collect all released material in a plastic lined container. Dispose off according to the local law and rules. Avoid leached substances to get into the earth, canalization or waters.	
<b>Section 7- Handling and Storage</b>	
<b>Handling</b>	The battery should not be opened, destroyed or incinerate, since they may leak or rupture and release to the environment the ingredients that they contain in the hermetically sealed container. Do not short circuit terminals, or over charge the battery, forced over-discharge, throw to fire. Do not crush or puncture the battery, or immerse in liquids.
<b>Storage</b>	Avoid mechanical or electrical abuse. Storage preferably in cool, dry and ventilated area, which is subject to little temperature change. Storage at high temperatures should be avoided. Do not place the battery near heating equipment, nor expose to direct sunlight for long periods.
<b>Other Precautions</b>	The battery may explode or cause burns, if disassembled, crushed or exposed to fire or high temperatures. Do not short or install with incorrect polarity.
<b>Section 8- Exposure Controls/Personal Protection</b>	
<b>Engineering Controls</b>	Use local exhaust ventilation or other engineering controls to control sources of dust, mist, fumes and vapor. Keep away from heat and open flame. Store in a cool, dry place.
<b>Personal Protective Equipment</b>	Respiratory Protection: Not necessary under normal conditions. Skin and body Protection: Not necessary under normal conditions, Wear suitable protective clothing and gloves if handling an open or leaking battery. Hand protection: Wear suitable gloves if handling an open or leaking battery. Eye Protection: Not necessary under normal conditions, Wear safety glasses if handling an open or leaking battery.
<b>Other Protective Equipment</b>	Have a safety shower and eye wash fountain readily available in the immediate work area.
<b>Hygiene Measures</b>	Do not eat, drink, or smoke in work area. Maintain good housekeeping.
<b>Section 9- Physical and Chemical Properties</b>	
<b>Form</b>	Solid
<b>Color</b>	Black
<b>Odour</b>	Not Applicable
<b>pH</b>	Not Applicable

## Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

<b>Melting point/freezing point</b>	Not Applicable
<b>Boiling Point and Boiling range</b>	Not Applicable
<b>Flash Point</b>	Not Applicable
<b>Upper/lower flammability or explosive limits</b>	Not Applicable
<b>Vapor Pressure</b>	Not Applicable
<b>Vapor Density</b>	Not Applicable
<b>Relative density</b>	Not Applicable
<b>Solubility in Water</b>	Not Applicable
<b>Auto-ignition temperature</b>	Not Applicable
<b>Decomposition temperature</b>	Not Applicable
<b>Evaporation rate</b>	Not Applicable
<b>Flammability (soil, gas)</b>	Not Applicable
<b>Viscosity</b>	Not Applicable
<b>Section 10- Stability and reactivity</b>	
<b>Stability</b>	The product is stable under conditions described Section 7
<b>Conditions to Avoid</b>	Heat above 70°C or incinerate. Deform, Mutilate, Crush, Disassemble, Overcharge, Short circuit, Expose over a long period to humid conditions.
<b>Incompatible Materials</b>	Oxidizing agents, acid, base.
<b>Hazardous Decomposition Products</b>	Carbon monoxide, carbon dioxide, lithium oxide fumes.
<b>Possibility of Hazardous Reaction</b>	Not Applicable
<b>Section 11 – Toxicological Information</b>	
<b>Irritation</b>	Risk of irritation occurs only if the cell is mechanically, thermally or electrically abused to the point of compromising the enclosure. If this occurs, irritation to the skin, eyes and respiratory tract may occur.
<b>Sensitization</b>	Not Applicable
<b>Neurological Effects</b>	Not Applicable
<b>Teratogenicity</b>	Not Applicable
<b>Reproductive Toxicity</b>	Not Applicable

## Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

<b>Mutagenicity (Genetic Effects)</b>	Not Applicable	
<b>Toxicologically Synergistic Materials</b>	Not Applicable	
<b>Section 12- Ecological Information</b>		
<b>Ecological Toxicity</b>	Not Applicable	
<b>Mobility in soil</b>	Not Applicable	
<b>Persistence and Degradability</b>	Not Applicable	
<b>Bioaccumulation potential</b>	Not Applicable	
<b>Other Adverse Effects</b>	Not Applicable	
<b>Section 13- Disposal Considerations</b>		
<b>Product disposal recommendation</b>	Observe local, state and federal laws and regulations.	
<b>Uncleaned packaging recommendation</b>	Disposal must be made according to official regulations	
<b>Section 14 – Transport Information</b>		
<b>Label for conveyance</b>	Lithium Battery Label	
<b>UN Number</b>	UN 3480 or UN 3171 or UN 3481	
<b>Transport hazard class(es)</b>	9	
<b>Packing group</b>	965	II
	952	---
	967	---
	966	II
<b>Marine pollutant</b>	No	
<b>UN Proper shipping name</b>	Lithium ion Batteries (Including lithium ion polymer batteries) Battery-powered vehicle or Battery-powered equipment Lithium ion Batteries packed with equipment (Including lithium ion polymer batteries) Lithium ion Batteries contained in equipments (Including lithium ion polymer batteries)	
<b>ICAO/IATA</b>	Can be shipped by air in accordance with international Civil Aviation Organization (ICAO), TI or International Air Transport Association (IATA) DGR 63 <sup>rd</sup> Packing Instructions Section IA of 965 or Section I of 966 967 or 952 appropriately.	

# Material Safety Data Sheet

REPORT NO.: KA2112161149AD100

<b>IMDG CODE</b>	International Maritime Dangerous Goods Code IMDG CODE (Amdt 40-20)
<b>ADR</b>	European Agreement concerning the International Carriage of Dangerous Goods by Road
<b>RID</b>	Regulations concerning the International Carriage of Dangerous Goods by Rail
The dangerous goods regulations require that each battery design be subject to tests contained in Section 38.3 of the UN Manual of Tests and Criteria prior to being offered for transport.	
<b>Section 15- Regulatory information</b>	
<p><b>Law information</b></p> <p>《Dangerous Goods Regulations》  《Recommendation on the Transport of Dangerous Goods Model Regulations》  《International Maritime Dangerous Goods》  《Technical Instructions for the Safe Transport of Dangerous Goods》  《Classification and code of dangerous Goods》  《Consumer Product Safety Act》 (CPSA)  《Federal Environmental Pollution Control Act》 (FEPCA)  《Resource Conservation and Recovery Act》 (RCRA)  《European Agreement concerning the International Carriage of Dangerous》  《Regulations concerning the International Carriage of Dangerous》</p> <p>In according with all Federal, State and local laws.</p>	
<b>Section 16- Other Information</b>	
<p>The information above is believed to be accurate and represents the best information currently available to us. However, Concorde makes no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. Although reasonable precautions have been taken in the preparation of the data contained herein, it is offered solely for your information, consideration and investigation. This material safety data sheet provides guidelines for the safe handling and use of this product; it does not and cannot advise on all possible situations, therefore, your specific use of this product should be evaluated to determine if additional precautions are required.</p>	

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## Lithium Battery UN38.3 Test Report

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ST/SG/AC.10/11/Rev.7/Subsection 38.3

**Sample name : Polymer Lithium-Ion Battery**

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Details information of the battery and the cell built in the battery, as following:

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